

## **Control Sacchrosomes without Cytochrome P450 Reductase**

### ***Product overview***

Catalogue Number CON001      Lot Number CON001-10-07

Membranes consist of a yeast microsomal fraction which is absent of any Cytochrome P450s or Cytochrome P450 Reductase enzyme.

### ***Product Data***

Pack Size	0.5 ml
Protein Concentration	10 mg.ml <sup>-1</sup>

### ***Protein Concentration***

Total protein was measured using a microplate-based Bradford assay method with BSA as a standard.

### ***Product Use***

For best stability thaw on ice, aliquot suitable quantities for your type of studies and store at -80°C.

Microsomes are supplied in a buffer containing Water, Tris, EDTA and Glycerol which are unlikely to interfere with most assays.

Studies indicate product stability at -80°C for at least 12 months.

### ***Safety***

This product is not suspected to contain any pathogenic or hazardous contents. However, since these properties have not been investigated handle with care in accordance to your normal laboratory practices.

This product is only intended for *in vitro* research use and is not licensed as a drug, therapy or diagnostic tool for humans or animals.

Sacchrosomes™ of DMBioscience Ltd, Patents pending, ©.